



BRAUN


BUA6150CEME

Quick reference guide

Please keep with unit / съхранявайте заедно с апарата / uchovávejte u přístroje / molimo čuvajte uz uređaj / kérjük, tartsa a termék mellett / przechowować razem z urządzeniem / الرجاء الاحتفاظ بها مع الوحدة / hranite skupaj z aparatom / Prosim, uchovajte spolu s prístrojom / Lütfen ürünle birlikte bulundurun



BUA6150CEME



GB **BG**



Irregular heart beat
Неправилен сърдечен ритъм



Pulse symbol
Символ за пулса

M:

Memory
Памет



Averaging function
Функция за осредняване



Low battery
Изтощена батерия

Error

Error (see owner's manual)
Грешка (вж. ръководството на потребителя)



Nepravidelný srdeční tep
Nepravilan rad srca



Symbol pulsu
Puls-Symbol



Paměť
Memorija



Výpočet průměru
Funkcija srednjih vrijednosti



Slabá baterie
Slaba baterija



Chyba (viz návod k obsluze)
Greška (pogledajte priručnik za korisnike)



Szabálytalan szívverés
Nieregularne bicie serca



Pulzus szimbóluma
Symbol tętna



Memória
Pamięć



Átlagérték-kiszámítási funkció
Funkcja uśredniania



Alacsony elemtöltés
Niski poziom na ładowania baterii



Hiba (lásd az útmutatót)
Błąd (patrz pod ręcznik użytkownika)



نبض قلب غير منتظم
Nepravilen srčni utrip



رمز النبض
Simbol za srčni utrip



الذاكرة
Pomnilnik



وظيفة المقاربة
Izračun povprečja



البطارية منخفضة
Baterija je skoraj prazna



خطأ (انظر دليل المالك)
Napaka (glejte uporabniški priročnik)



Nepravidelná srdcová frekvencia
Düzensiz nabız



Symbol pulzu
Nabız simgesi



Pamäť
Hafıza



Funkcia spriemerovania
Ortalama alma işlevi



Slabá batéria
Zayıf pil



Chyba (pozri príručku majiteľa)
Hata (Bknz. kullanma kılavuzu)

This blood pressure monitor is manufactured under licence to the 'Braun' trademark.

'Braun' is a registered trademark of Braun GmbH, Kronberg, Germany.

Helen of Troy

CREATIVE DEPARTMENT ARTWORK SPECIFICATIONS

Brand: Braun
Category: BPM

Model: **BUA6150CEME**
Artwork Part #: A002608R0
Die Line Part #: N/A
Subject: Quick Reference Guide
Region: EMEA
Flat Size: MM: W58 x H225
Folded Size: MM: W58 x H45

Scale: 1/1
Material: 50# matte text

Page Count: 2
Revision: 3
Date: **24APR18**
Release Date: 14MAY18
Re-release Date:

Colors:

Dielines (Do not print)

Spot Colors

Cyan
0%

Magenta
0%

Yellow
0%

Black
100%



Special Instructions:

**Quality Requirement of Artwork and Quality Clarification Process of Artwork Printing:
Meet Eng-QS-06&02**

Kaz USA, Inc., a Helen of Troy Company
Creative Services
Marlborough, MA 01752, USA
+1 508 490 7000